



PCSK7 Antibody, Biotin conjugated

Product Code	CSB-PA017646LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16549
Immunogen	Recombinant Human Proprotein convertase subtilisin/kexin type 7 protein (189-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proprotein convertase subtilisin/kexin type 7 (EC 3.4.21.-) (Lymphoma proprotein convertase) (Prohormone convertase 7) (Proprotein convertase 7) (PC7) (Proprotein convertase 8) (PC8) (hPC8) (Subtilisin/kexin-like protease PC7), PCSK7, LPC PC7 PC8 SPC7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCSK7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.