



PLA1A Antibody, Biotin conjugated

Product Code	CSB-PA018084LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53H76
Immunogen	Recombinant Human Phospholipase A1 member A protein (293-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase A1 member A (EC 3.1.1.-) (Phosphatidylserine-specific phospholipase A1) (PS-PLA1), PLA1A, NMD PSPLA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLA1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.