



POLD2 Antibody, Biotin conjugated

Product Code	CSB-PA018304LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49005
Immunogen	Recombinant Human DNA polymerase delta subunit 2 protein (245-430AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase delta subunit 2 (DNA polymerase delta subunit p50), POLD2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.