




POLG Antibody

Product Code	CSB-PA018313LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54098
Immunogen	Recombinant Human DNA polymerase subunit gamma-1 protein (446-590AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase subunit gamma-1 (EC 2.7.7.7) (Mitochondrial DNA polymerase catalytic subunit) (PolG-alpha), POLG, MDP1 POLG1 POLGA
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLG

Image



Western Blot

Positive WB detected in: 293T whole cell lysate

All lanes: POLG antibody at 4.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 140 kDa

Observed band size: 140 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.