



PTPN3 Antibody, Biotin conjugated

Product Code	CSB-PA019040LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26045
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 3 protein (365-505AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 3 (EC 3.1.3.48) (Protein-tyrosine phosphatase H1) (PTP-H1), PTPN3, PTPH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.