



RFX1 Antibody, Biotin conjugated

Product Code	CSB-PA019610LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22670
Immunogen	Recombinant Human MHC class II regulatory factor RFX1 protein (36-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MHC class II regulatory factor RFX1 (Enhancer factor C) (EF-C) (Regulatory factor X 1) (RFX) (Transcription factor RFX1), RFX1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RFX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.