



RFXANK Antibody, Biotin conjugated

Product Code	CSB-PA019617LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14593
Immunogen	Recombinant Human DNA-binding protein RFXANK protein (3-84AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein RFXANK (Ankyrin repeat family A protein 1) (Regulatory factor X subunit B) (RFX-B) (Regulatory factor X-associated ankyrin-containing protein), RFXANK, ANKRA1 RFXB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	RFXANK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.