



RPS29 Antibody, FITC conjugated

Product Code	CSB-PA020432LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62273
Immunogen	Recombinant Human 40S ribosomal protein S29 protein (2-56AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S29 (Small ribosomal subunit protein uS14), RPS29
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.