



ENTR1 Antibody, Biotin conjugated

Product Code	CSB-PA020897LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C92
Immunogen	Recombinant Human Endosome-associated-trafficking regulator 1 protein (121-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endosome-associated-trafficking regulator 1 (Antigen NY-CO-3) (Serologically defined colon cancer antigen 3), ENTR1, SDCCAG3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ENTR1