



ENTR1 Antibody, Biotin conjugated

Product Code	CSB-PA020897LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C92
Immunogen	Recombinant Human Endosome-associated-trafficking regulator 1 protein (121-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endosome-associated-trafficking regulator 1 (Antigen NY-CO-3) (Serologically defined colon cancer antigen 3), ENTR1, SDCCAG3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ENTR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.