




SFRP5 Antibody

Product Code	CSB-PA021141LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T4F7
Immunogen	Recombinant Human Secreted frizzled-related protein 5 protein (169-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 5 (sFRP-5) (Frizzled-related protein 1b) (FRP-1b) (Secreted apoptosis-related protein 3) (SARP-3), SFRP5, FRP1B SARP3
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SFRP5

Image



IHC image of CSB-PA021141LA01HU diluted at 1:200 and staining in paraffin-embedded human cervical cancer performed on a Leica Bond™ system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.