



SFRP5 Antibody, Biotin conjugated

Product Code	CSB-PA021141LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T4F7
Immunogen	Recombinant Human Secreted frizzled-related protein 5 protein (169-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 5 (sFRP-5) (Frizzled-related protein 1b) (FRP-1b) (Secreted apoptosis-related protein 3) (SARP-3), SFRP5, FRP1B, SARP3
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SFRP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.