



SLAMF6 Antibody, FITC conjugated

Product Code	CSB-PA021368LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DU3
Immunogen	Recombinant Human SLAM family member 6 protein (248-331AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLAM family member 6 (Activating NK receptor) (NK-T-B-antigen) (NTB-A) (CD antigen CD352), SLAMF6, KALI
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLAMF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.