



SLC19A2 Antibody, FITC conjugated

Product Code	CSB-PA021430LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60779
Immunogen	Recombinant Human Thiamine transporter 1 protein (209-285AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiamine transporter 1 (ThTr-1) (ThTr1) (Solute carrier family 19 member 2) (Thiamine carrier 1) (TC1), SLC19A2, THT1 TRMA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC19A2