



SLC19A2 Antibody, Biotin conjugated

Product Code	CSB-PA021430LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60779
Immunogen	Recombinant Human Thiamine transporter 1 protein (209-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thiamine transporter 1 (ThTr-1) (ThTr1) (Solute carrier family 19 member 2) (Thiamine carrier 1) (TC1), SLC19A2, THT1 TRMA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC19A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.