



SLC26A4 Antibody, HRP conjugated

Product Code	CSB-PA021527LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43511
Immunogen	Recombinant Human Pendrin protein (592-655AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pendrin (Sodium-independent chloride/iodide transporter) (Solute carrier family 26 member 4), SLC26A4, PDS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC26A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.