



SLC2A4RG Antibody, Biotin conjugated

Product Code	CSB-PA021558LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR83
Immunogen	Recombinant Human SLC2A4 regulator protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLC2A4 regulator (GLUT4 enhancer factor) (GEF) (Huntington disease gene regulatory region-binding protein 1) (HDBP-1), SLC2A4RG, HDBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC2A4RG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.