



SLC7A1 Antibody, FITC conjugated

Product Code	CSB-PA021709LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30825
Immunogen	Recombinant Human High affinity cationic amino acid transporter 1 protein (430-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity cationic amino acid transporter 1 (CAT-1) (CAT1) (Ecotropic retroviral leukemia receptor homolog) (Ecotropic retrovirus receptor homolog) (Solute carrier family 7 member 1) (System Y+ basic amino acid transporter), SLC7A1, ATRC1 ERR REC1L
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC7A1