



SLC7A1 Antibody, FITC conjugated

Product Code	CSB-PA021709LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30825
Immunogen	Recombinant Human High affinity cationic amino acid transporter 1 protein (430-492AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity cationic amino acid transporter 1 (CAT-1) (CAT1) (Ecotropic retroviral leukemia receptor homolog) (Ecotropic retrovirus receptor homolog) (Solute carrier family 7 member 1) (System Y+ basic amino acid transporter), SLC7A1, ATRC1 ERR REC1L
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC7A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.