




SOX8 Antibody

Product Code	CSB-PA022436LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57073
Immunogen	Recombinant Human Transcription factor SOX-8 protein (292-386AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor SOX-8, SOX8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SOX8

Image



Western Blot

Positive WB detected in: U87 whole cell lysate
 All lanes: SOX8 antibody at 5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 48 kDa
 Observed band size: 48 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.