



SPOPL Antibody, FITC conjugated

Product Code	CSB-PA022602LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IQ16
Immunogen	Recombinant Human Speckle-type POZ protein-like protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Speckle-type POZ protein-like (HIB homolog 2) (Roadkill homolog 2), SPOPL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPOPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.