



SYNJ1 Antibody, Biotin conjugated

Product Code	CSB-PA023016LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43426
Immunogen	Recombinant Human Synaptojanin-1 protein (1030-1261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptojanin-1 (EC 3.1.3.36) (Synaptic inositol 1,4,5-trisphosphate 5-phosphatase 1), SYNJ1, KIAA0910
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SYNJ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.