



TFB1M Antibody, Biotin conjugated

Product Code	CSB-PA023420LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVM0
Immunogen	Recombinant Human Dimethyladenosine transferase 1, mitochondrial protein (59-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dimethyladenosine transferase 1, mitochondrial (EC 2.1.1.-) (Mitochondrial 12S rRNA dimethylase 1) (Mitochondrial transcription factor B1) (h-mtTFB) (h-mtTFB1) (hTFB1M) (mtTFB1) (S-adenosylmethionine-6-N', N'-adenosyl(rRNA) dimethyltransferase 1), TFB1M
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TFB1M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.