




TMSB4X Antibody

Product Code	CSB-PA023934LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62328
Immunogen	Recombinant Human Thymosin beta-4 protein (6-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymosin beta-4 (T beta-4) (Fx) [Cleaved into: Hematopoietic system regulatory peptide (Seraspenide)], TMSB4X, TB4X THYB4 TMSB4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMSB4X

Image



Immunohistochemistry of paraffin-embedded human lymph node tissue using CSB-PA023934LA01HU at dilution of 1:100