



TRAPPC6A Antibody, Biotin conjugated

Product Code	CSB-PA024232LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75865
Immunogen	Recombinant Human Trafficking protein particle complex subunit 6A protein (1-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 6A (TRAPP complex subunit 6A), TRAPPC6A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRAPPC6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.