



TRAPPC6A Antibody, Biotin conjugated

Product Code	CSB-PA024232LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75865
Immunogen	Recombinant Human Trafficking protein particle complex subunit 6A protein (1-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 6A (TRAPP complex subunit 6A), TRAPPC6A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRAPPC6A