



UNC5B Antibody, Biotin conjugated

Product Code	CSB-PA025629LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZJ1
Immunogen	Recombinant Human Netrin receptor UNC5B protein (404-517AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Netrin receptor UNC5B (Protein unc-5 homolog 2) (Protein unc-5 homolog B) (p53-regulated receptor for death and life protein 1) (p53RDL1), UNC5B, P53RDL1 UNC5H2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC5B