



UNC5B Antibody, Biotin conjugated

Product Code	CSB-PA025629LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZJ1
Immunogen	Recombinant Human Netrin receptor UNC5B protein (404-517AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Netrin receptor UNC5B (Protein unc-5 homolog 2) (Protein unc-5 homolog B) (p53-regulated receptor for death and life protein 1) (p53RDL1), UNC5B, P53RDL1 UNC5H2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC5B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.