



WNT9B Antibody, HRP conjugated

Product Code	CSB-PA026146LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14905
Immunogen	Recombinant Human Protein Wnt-9b protein (213-324AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wnt-9b (Protein Wnt-14b) (Protein Wnt-15), WNT9B, WNT14B WNT15
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	WNT9B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.