



XIRP2 Antibody, FITC conjugated

Product Code	CSB-PA026196LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4UGR9
Immunogen	Recombinant Human Xin actin-binding repeat-containing protein 2 protein (2861-3004AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Xin actin-binding repeat-containing protein 2 (Beta-xin) (Cardiomyopathy-associated protein 3) (Xeplin), XIRP2, CMYA3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	XIRP2