



XIRP2 Antibody, Biotin conjugated

Product Code	CSB-PA026196LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4UGR9
Immunogen	Recombinant Human Xin actin-binding repeat-containing protein 2 protein (2861-3004AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Xin actin-binding repeat-containing protein 2 (Beta-xin) (Cardiomyopathy-associated protein 3) (Xeplin), XIRP2, CMYA3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	XIRP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.