



ZNF703 Antibody, Biotin conjugated

Product Code	CSB-PA026966LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7S9
Immunogen	Recombinant Human Zinc finger protein 703 protein (218-286AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 703 (Zinc finger elbow-related proline domain protein 1), ZNF703, ZEPP01 ZPO1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF703
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.