



# ADAM8 Antibody, Biotin conjugated

Product Code	CSB-PA10259D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78325
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 8 protein (677-824AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 8 (ADAM 8) (EC 3.4.24.-) (Cell surface antigen MS2) (CD antigen CD156a), ADAM8, MS2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ADAM8