



ADAM8 Antibody, Biotin conjugated

Product Code	CSB-PA10259D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78325
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 8 protein (677-824AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 8 (ADAM 8) (EC 3.4.24.-) (Cell surface antigen MS2) (CD antigen CD156a), ADAM8, MS2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ADAM8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.