



MGAT3 Antibody, HRP conjugated

Product Code	CSB-PA22969B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09327
Immunogen	Recombinant Human Beta-1,4-mannosyl-glycoprotein 4-beta-N-acetylglucosaminyltransferase protein (435-533AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-1,4-mannosyl-glycoprotein 4-beta-N-acetylglucosaminyltransferase (EC 2.4.1.144) (N-glycosyl-oligosaccharide-glycoprotein N-acetylglucosaminyltransferase III) (GNT-III) (GlcNAc-T III) (N-acetylglucosaminyltransferase III), MGAT3, GGNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MGAT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.