



UBIAD1 Antibody, Biotin conjugated

Product Code	CSB-PA23029D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Z9
Immunogen	Recombinant Human UbiA prenyltransferase domain-containing protein 1 protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UbiA prenyltransferase domain-containing protein 1 (EC 2.5.1.-) (Transitional epithelial response protein 1), UBIAD1, TERE1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	UBIAD1