



IGSF9 Antibody, HRP conjugated

| | |
|----------------------------|-------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA23039B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9P2J2 |
| Immunogen | Recombinant Human Protein turtle homolog A protein (132-298AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein turtle homolog A (Immunoglobulin superfamily member 9A) (IgSF9A), IGSF9, IGSF9A KIAA1355 NRT1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | IGSF9 |