



IGSF9 Antibody, Biotin conjugated

Product Code	CSB-PA23039D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2J2
Immunogen	Recombinant Human Protein turtle homolog A protein (132-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein turtle homolog A (Immunoglobulin superfamily member 9A) (IgSF9A), IGSF9, IGSF9A KIAA1355 NRT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	IGSF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.