



DPY19L1 Antibody, FITC conjugated

Product Code	CSB-PA23079C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2PZI1
Immunogen	Recombinant Human Probable C-mannosyltransferase DPY19L1 protein (513-665AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable C-mannosyltransferase DPY19L1 (EC 2.4.1.-) (Dpy-19-like protein 1) (Protein dpy-19 homolog 1), DPY19L1, GA0500 KIAA0877
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPY19L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.