



# DPY19L1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA23079D0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q2PZI1
<b>Immunogen</b>	Recombinant Human Probable C-mannosyltransferase DPY19L1 protein (513-665AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Probable C-mannosyltransferase DPY19L1 (EC 2.4.1.-) (Dpy-19-like protein 1) (Protein dpy-19 homolog 1), DPY19L1, GA0500 KIAA0877
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	DPY19L1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.