



Alpha-insect toxin LqhaIT Antibody

Product Code	CSB-PA322958LA01LDS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17728
Immunogen	Recombinant Leiurus quinquestriatus hebraeus Alpha-insect toxin LqhaIT protein (20-85AA)
Raised In	Rabbit
Species Reactivity	Leiurus quinquestriatus hebraeus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-insect toxin LqhaIT (Lqh-alpha-IT) (Alpha-IT)
Immunogen Species	Leiurus quinquestriatus hebraeus (Yellow scorpion)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.