



CEL Antibody, Biotin conjugated

Product Code	CSB-PA325708LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19835
Immunogen	Recombinant Human Bile salt-activated lipase protein (661-753AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bile salt-activated lipase (BAL) (EC 3.1.1.13) (EC 3.1.1.3) (Bile salt-stimulated lipase) (BSSL) (Bucelipase) (Carboxyl ester lipase) (Cholesterol esterase) (Pancreatic lysophospholipase) (Sterol esterase), CEL, BAL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.