




BZLF2 Antibody

Product Code	CSB-PA365876LA01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03205
Immunogen	Recombinant Epstein-Barr virus Glycoprotein 42 protein (34-223AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycoprotein 42 (gp42) [Cleaved into: Soluble gp42]
Immunogen Species	Epstein-Barr virus
Research Area	Others
Target Names	BZLF2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: BZLF2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 38 kDa

Observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.