



BZLF2 Antibody, HRP conjugated

Product Code	CSB-PA365876LB01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03205
Immunogen	Recombinant Epstein-Barr virus Glycoprotein 42 protein (34-223AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycoprotein 42 (gp42) [Cleaved into: Soluble gp42]
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Research Area	Others
Target Names	BZLF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.