



RBFOX3 Antibody, Biotin conjugated

Product Code	CSB-PA407385LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NFN3
Immunogen	Recombinant Human RNA binding protein fox-1 homolog 3 protein (1-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA binding protein fox-1 homolog 3 (Fox-1 homolog C) (Neuronal nuclei antigen) (NeuN antigen), RBFOX3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RBFOX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.