



ST3GAL4 Antibody, Biotin conjugated

Product Code	CSB-PA612121LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q11206
Immunogen	Recombinant Human CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 4 protein (39-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 4 (Alpha 2,3-ST 4) (Beta-galactoside alpha-2,3-sialyltransferase 4) (EC 2.4.99.2) (EC 2.4.99.4) (Alpha 2,3-sialyltransferase IV) (Gal-NAc6S) (Gal-beta-1,4-GalNAc-alpha-2,3-sialyltransferase) (SAT-3) (ST-4) (ST3Gal IV) (ST3GalIV) (ST3GalA.2) (STZ) (Sialyltransferase 4C) (SIAT4-C), ST3GAL4, CGS23 NANTA3 SIAT4C STZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ST3GAL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.