



PTP4A2 Antibody, Biotin conjugated

Product Code	CSB-PA613277LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12974
Immunogen	Recombinant Human Protein tyrosine phosphatase type IVA 2 protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein tyrosine phosphatase type IVA 2 (EC 3.1.3.48) (HU-PP-1) (OV-1) (PTP(CAAXII)) (Protein-tyrosine phosphatase 4a2) (Protein-tyrosine phosphatase of regenerating liver 2) (PRL-2), PTP4A2, PRL2 PTPCAAX2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTP4A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.