



DSC3 Antibody, HRP conjugated

Product Code	CSB-PA613507LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14574
Immunogen	Recombinant Human Desmocollin-3 protein (712-896AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Desmocollin-3 (Cadherin family member 3) (Desmocollin-4) (HT-CP), DSC3, CDHF3 DSC4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DSC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.