



CKAP5 Antibody, HRP conjugated

Product Code	CSB-PA617899LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14008
Immunogen	Recombinant Human Cytoskeleton-associated protein 5 protein (1802-2019AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoskeleton-associated protein 5 (Colonic and hepatic tumor overexpressed gene protein) (Ch-TOG), CKAP5, KIAA0097
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CKAP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.