



GINS1 Antibody, FITC conjugated

Product Code	CSB-PA617926LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14691
Immunogen	Recombinant Human DNA replication complex GINS protein PSF1 protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA replication complex GINS protein PSF1 (GINS complex subunit 1), GINS1, KIAA0186 PSF1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GINS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.